



6 Dry column	Centrifuge the TG Mini Column at full speed for an additional 3 min to dry the TG Mini Column.
7 DNA Elution	<b>1. Add 100 µl of preheated Elution Buffer or ddH2O (pH 7.5-9.0)</b> to the membrane of the TG Mini Column. Stand the TG Mini Column for 3 min. 2. Centrifuge at full speed for 2 min to elute DNA.

Special Protocol: The sample preparation For Animal Cultured Cells

Additional requirement	1. RNase A (optional) 2. 96~100% ethanol 3. trypsin or cell scraper (for monolayer cell ) 4. PBS
Method	<b>Harvest cells</b> <b>1. For Cells grown in suspension</b> a. Transfer the appropriate number of cell ( up to 1 x 10 <sup>7</sup> ) to a microcentrifuge tube. b. Centrifuge at 300 x g for 5 min. Discard supernatant carefully and completely.  <b>For Cells grown in monolayer</b> a. Detach cells from the dish or flask by trypsinization or using a cell scraper. Transfer the appropriate number of cell ( up to 1 x 10 <sup>7</sup> ) to a microcentrifuge tube. b. Centrifuge at 300 x g for 5 min. Discard supernatant carefully and completely.  2. Resuspend cell pellet in PBS to a final volume of <b>200 µl</b> . 3. Follow the Animal Tissuel Protocol starting from step 2.1 Lysis.

Special Protocol: The sample preparation For Blood

Additional requirement	1. RNase A (optional) 2. 96~100% ethanol 3. PBS
Method	1. Transfer up to <b>200 µl</b> sample ( whole blood, serum, plasma, body fluids, buffy coat) to a microcentrifuge tube. • <b>If the sample volume is less than 200 µl , add the appropriate volume of PBS.</b> <b>Optional:</b> If RNA-free genomic DNA is required, <b>add 4 µl of 100 mg/mL RNase A (not provided).</b> Mix thoroughly by vortexing and incubate at room temperature for 2 min.  <b>2. Add 20 µl Proteinase K</b> to the sample, and then <b>add 200 µl TG2 Buffer</b> to the sample. Mix thoroughly by pulse-vortexing. Incubate at 60°C for 30 min. Vrotex occasionally during incubation. 3. Incubate at 70°C for 10 min. 4. Follow the Animal Tissuel Protocol starting from step 3. Ethanol Dilution.

Special Protocol: The sample preparation For Bacteria

Additional requirement	1. RNase A (optional) 2. 96~100% ethanol 3. For Gram-positive bacteria: lysozyme reaction solution (20 mg/mL lysozyme; 20 mM Tris-HCl, pH 8.0; 2mM EDTA; 1.2 % Triton)
Method	<b>1. For bacterial cultures</b> a. Transfer 1 mL well-grown bacterial culture to a microcentrifuge tube (not provided). b. Descend the cells by centrifuging at full speed for 2 min and discard supernatant completely. c. Follow the Animal Tissue Protocol starting from step 2.1 Lysis.  <b>For bacterial in biological fluids</b> a. Collect cells by centrifuging biological fluids at 7,500 rpm (5,000 x g) for 10 min and discard supernatant completely. b. Follow the Animal Tissue Protocol starting from step 2.1 Lysis.  <b>For bacteria from eye, nasal, pharyngeal, or other swabs</b> a. Soak the swabs in 2 mL PBS at room temperature for 2- 3 hr. b. Collect cells by centrifuging at 7,500 rpm (5,000 x g) for 10 min and discard supernatant completely. c. Follow the Animal Tissue Protocol starting from step 2.1 Lysis.  <b>For Gram-positive bacteria</b> <b>HINT: Set dry or water baths: one to 37°C , another to 60°C and the other to 95°C .</b> a. Transfer 1 mL well-grown bacterial culture to a microcentrifuge tube (not provided). b. Descend the cells by centrifuging at full speed for 2 min and discard supernatant completely. c. Resuspend the cell pellet in <b>200 µl lysozyme reaction solution</b> (20 mg/mL lysozyme; 20 mM Tris-HCl, pH 8.0; 2mM EDTA; 1.2 % Triton). Incubate at 37°C for 30~60 min.  <b>Optional:</b> If RNA-free genomic DNA is required, add 4 µl of 100 mg/mL RNase A (not provided). Mix thoroughly by vortexing and incubate at room temperature for 2 min.  <b>2. Add 20 µl Proteinase K</b> to the sample, and then <b>add 200 µl TG2 Buffer</b> to the sample. Mix thoroughly by pulse-vortexing. Incubate at 60°C for 30 min and vrotex occasionally during incubation.  <b>Optional:</b> If required, incubate at 95°C for 15 min to inactivate pathogens. <b>Note:</b> that incubation at 95°C can lead to some DNA degradation.  3. Follow the Animal Tissue Protocol starting from step 3. Ethanol Dilution.

Special Protocol: The sample preparation For Yeast

Additional requirement	1. xylene & ethanol (96~100%) 2. liquid nitrogen 3. Mortar a rotor-stator homogenizer or a 20-G needle syringe 4. β-Mercaptoethanol 5. 70% RNase-free ethanol
Method	1. Transfer 3 mL log-phase (OD600 = 10) yeast culture to a microcentrifuge tube (not provided). 2. Descend the cells by centrifuging at 7,500 rpm (5,000 x g) for 10 min. Discard supernatant completely. 3. Resuspend the cell pellet in <b>600 µl sorbitol buffer</b> (1M sorbitol; 100 mM EDTA; 14 mM β-mercaptoethanol). <b>Add 200 U zymolase or lyticase</b> and incubate at 30°C for 30 min. 4. Centrifuge at 7,500 rpm (5,000 x g) for 5 min. Remove supernatant by pipetting. 5. Follow the Animal Tissue Protocol starting from step 2.1 Lysis.

Special Protocol: The sample preparation For Dried Blood Spot

Additional equipment	1. RNase A (optional) 2. 96~100% ethanol
Method	1. Cut the filter paper (e.g. S&S903) with dried blood spot into a microcentrifuge tube. Add <b>200 µl TG1 Buffer</b> and incubate at 85°C for 10 min. <b>2. Add 20 µl Proteinase K</b> to the sample mixture. Mix thoroughly by vortexing. Incubate at 60°C for 1 hr. Vrotex occasionally during incubation. 3. Follow the Animal Tissue Protocol starting from step 2.2 Lysis.

Special Protocol: The sample preparation For Fixed Tissue

Additional equipment	1. RNase A (optional) 2. 96~100% ethanol 3. Xylene - for paraffin-embedded tissues
Method	<b>For paraffin-embedded tissues</b> 1. Cut up to <b>25 mg</b> paraffin-embedded tissue sample to a microcentrifuge tube (not provided). <b>2. Add 1 mL xylene</b> , mix well and incubate at room temperature for 30 min. 3. Centrifuge at full speed for 5 min. Remove supernatant by pipetting. <b>4. Add 1 mL ethanol (96- 100 %)</b> to the deparaffined tissue, mix gently by vortexing. 5. Centrifuge at full speed for 3 min. Remove supernatant by pipetting. 6. Repeat step 4 and 5. 7. Incubate at 37°C for 10 ~15 min to evaporate ethanol residue completely. 8. Grind the tissue sample by micropestle or liquid nitrogen and follow the Animal Tissue Protocol starting from step 2.1 Lysis.  <b>For formalin-fixed tissues</b> 1. Wash <b>25 mg</b> tissue sample twice with <b>1 mL</b> PBS to remove formalin. 2. Grind the tissue sample by micropestle or liquid nitrogen and follow the Animal Tissue Protocol starting from step 2.1 Lysis.

TEL : 886-2-2228-8950  
FAX : 886-2-2228-8955  
Mail : [info@lis-tek.com](mailto:info@lis-tek.com)  
4F., No.926, Zhongzheng Rd., Zhonghe  
Dist., New Taipei City 235, Taiwan (R.O.C.)

[www.lis-bio.com/zh-TW](http://www.lis-bio.com/zh-TW)

